

```
import pandas as pd
import scipy.signal as scisig
import os
import numpy as np
```

```
def get_user_input(prompt):
    try:
        return raw_input(prompt)
    except NameError:
        return input(prompt)
```

```
def getInputLoadFile():
    """Asks user for type of file and file path. Loads corresponding data.
```

OUTPUT:

```
    data: DataFrame, index is a list of timestamps at 8Hz, columns include
           AccelZ, AccelY, AccelX, Temp, EDA, filtered_eda
    """
```

```
print("Please enter information about your EDA file... ")
dataType = get_user_input("\tData Type (e4, q, shimmer, or misc): ")
if dataType=='q':
    filepath = get_user_input("\tFile path: ")
    filepath_confirm = filepath
    data = loadData_Qsensor(filepath)
elif dataType=='e4':
    filepath = get_user_input("\tPath to E4 directory: ")
    filepath_confirm = os.path.join(filepath, "EDA.csv")
```

```

    data = loadData_E4(filepath)
elif dataType=='shimmer':
    filepath = get_user_input("\tFile path: ")
    filepath_confirm = filepath
    data = loadData_shimmer(filepath)
elif dataType=="misc":
    filepath = get_user_input("\tFile path: ")
    filepath_confirm = filepath
    data = loadData_misc(filepath)
else:
    print("Error: not a valid file choice")

return data, filepath_confirm

def getOutputPath():
    print("")
    print("Where would you like to save the computed output file?")
    outfile = get_user_input('\tFile name: ')
    outputPath = get_user_input("\tFile directory (./ for this directory): ")
    fullOutputPath = os.path.join(outputPath,outfile)
    if fullOutputPath[-4:] != '.csv':
        fullOutputPath = fullOutputPath+'.csv'
    return fullOutputPath

def loadData_Qsensor(filepath):
    """
    This function loads the Q sensor data, uses a lowpass butterworth filter on the EDA signal

    Note: currently assumes sampling rate of 8hz, 16hz, 32hz; if sampling rate is 16hz or 32hz the signal is
    downsampled

```

INPUT:

filepath: string, path to input file

OUTPUT:

data: DataFrame, index is a list of timestamps at 8Hz, columns include AccelZ, AccelY, AccelX, Temp, EDA, filtered\_eda

'''

# Get header info

try:

header\_info = pd.io.parsers.read\_csv(filepath, nrows=5)

except IOError:

print("Error!! Couldn't load file, make sure the filepath is correct and you are using a csv from the q sensor software\n\n")

return

# Get sample rate

sampleRate = int((header\_info.iloc[3,0]).split(":")[1].strip())

# Get the raw data

data = pd.io.parsers.read\_csv(filepath, skiprows=7)

data = data.reset\_index()

# Reset the index to be a time and reset the column headers

data.columns = ['AccelZ', 'AccelY', 'AccelX', 'Battery', 'Temp', 'EDA']

# Get Start Time

startTime = pd.to\_datetime(header\_info.iloc[4,0][12:-10])

```

# Make sure data has a sample rate of 8Hz
data = interpolateDataTo8Hz(data,sampleRate,startTime)

# Remove Battery Column
data = data[['AccelZ','AccelY','AccelX','Temp','EDA']]

# Get the filtered data using a low-pass butterworth filter (cutoff:1hz, fs:8hz, order:6)
data['filtered_eda'] = butter_lowpass_filter(data['EDA'], 1.0, 8, 6)

return data

def _loadSingleFile_E4(filepath,list_of_columns, expected_sample_rate,freq):
    # Load data
    data = pd.read_csv(filepath)

    # Get the startTime and sample rate
    startTime = pd.to_datetime(float(data.columns.values[0]),unit="s")
    sampleRate = float(data.iloc[0][0])
    data = data[data.index!=0]
    data.index = data.index-1

    # Reset the data frame assuming expected_sample_rate
    data.columns = list_of_columns
    if sampleRate != expected_sample_rate:
        print('ERROR, NOT SAMPLED AT {0}HZ. PROBLEMS WILL OCCUR\n'.format(expected_sample_rate))

    # Make sure data has a sample rate of 8Hz
    data = interpolateDataTo8Hz(data,sampleRate,startTime)

```

```
return data
```

```
def loadData_E4(filepath):
```

```
    # Load EDA data
```

```
    eda_data = _loadSingleFile_E4(os.path.join(filepath, 'EDA.csv'), ["EDA"], 4, "250L")
```

```
    # Get the filtered data using a low-pass butterworth filter (cutoff:1hz, fs:8hz, order:6)
```

```
    eda_data['filtered_eda'] = butter_lowpass_filter(eda_data['EDA'], 1.0, 8, 6)
```

```
    # Load ACC data
```

```
    acc_data =
```

```
    _loadSingleFile_E4(os.path.join(filepath, 'ACC.csv'), ["AccelX", "AccelY", "AccelZ"], 32, "31250U")
```

```
    # Scale the accelerometer to +-2g
```

```
    acc_data[["AccelX", "AccelY", "AccelZ"]] = acc_data[["AccelX", "AccelY", "AccelZ"]]/64.0
```

```
    # Load Temperature data
```

```
    temperature_data = _loadSingleFile_E4(os.path.join(filepath, 'TEMP.csv'), ["Temp"], 4, "250L")
```

```
    data = eda_data.join(acc_data, how='outer')
```

```
    data = data.join(temperature_data, how='outer')
```

```
    # E4 sometimes records different length files - adjust as necessary
```

```
    min_length = min(len(acc_data), len(eda_data), len(temperature_data))
```

```
    return data[:min_length]
```

```
def loadData_shimmer(filepath):
```

```
    data = pd.read_csv(filepath, sep='\t', skiprows=(0,1))
```

```

orig_cols = data.columns
rename_cols = {}

for search, new_col in [['Timestamp', 'Timestamp'],
                        ['Accel_LN_X', 'AccelX'], ['Accel_LN_Y', 'AccelY'], ['Accel_LN_Z', 'AccelZ'],
                        ['Skin_Conductance', 'EDA']]:
    orig = [c for c in orig_cols if search in c]
    if len(orig) == 0:
        continue
    rename_cols[orig[0]] = new_col

data.rename(columns=rename_cols, inplace=True)

# TODO: Assuming no temperature is recorded
data['Temp'] = 0

# Drop the units row and unnecessary columns
data = data[data['Timestamp'] != 'ms']
data.index = pd.to_datetime(data['Timestamp'], unit='ms')
data = data[['AccelZ', 'AccelY', 'AccelX', 'Temp', 'EDA']]

for c in ['AccelZ', 'AccelY', 'AccelX', 'Temp', 'EDA']:
    data[c] = pd.to_numeric(data[c])

# Convert to 8Hz
data = data.resample("125L").mean()
data.interpolate(inplace=True)

# Get the filtered data using a low-pass butterworth filter (cutoff:1hz, fs:8hz, order:6)

```

```
data['filtered_eda'] = butter_lowpass_filter(data['EDA'], 1.0, 8, 6)
```

```
return data
```

```
def loadData_getColNames(data_columns):
```

```
    print("Here are the data columns of your file: ")
```

```
    print(data_columns)
```

```
    # Find the column names for each of the 5 data streams
```

```
    colnames = ['EDA data', 'Temperature data', 'Acceleration X', 'Acceleration Y', 'Acceleration Z']
```

```
    new_colnames = [",", ",", ",", ",", ","]
```

```
    for i in range(len(new_colnames)):
```

```
        new_colnames[i] = get_user_input("Column name that contains "+colnames[i]+": ")
```

```
        while (new_colnames[i] not in data_columns):
```

```
            print("Column not found. Please try again")
```

```
            print("Here are the data columns of your file: ")
```

```
            print(data_columns)
```

```
        new_colnames[i] = get_user_input("Column name that contains "+colnames[i]+": ")
```

```
    # Get user input on sample rate
```

```
    sampleRate = get_user_input("Enter sample rate (must be an integer power of 2): ")
```

```
    while (sampleRate.isdigit()==False) or (np.log(int(sampleRate))/np.log(2) !=  
np.floor(np.log(int(sampleRate))/np.log(2))):
```

```
        print("Not an integer power of two")
```

```
        sampleRate = get_user_input("Enter sample rate (must be a integer power of 2): ")
```

```
    sampleRate = int(sampleRate)
```

```

# Get user input on start time

startTime = pd.to_datetime(get_user_input("Enter a start time (format: YYYY-MM-DD HH:MM:SS): "))

while type(startTime)!=str:

    print("Not a valid date/time")

    startTime = pd.to_datetime(get_user_input("Enter a start time (format: YYYY-MM-DD HH:MM:SS):
"))

return sampleRate, startTime, new_colnames


def loadData_misc(filepath):

    # Load data

    data = pd.read_csv(filepath)


    # Get the correct colnames

    sampleRate, startTime, new_colnames = loadData_getColNames(data.columns.values)


    data.rename(columns=dict(zip(new_colnames,['EDA','Temp','AccelX','AccelY','AccelZ'])), inplace=True)

    data = data[['AccelZ','AccelY','AccelX','Temp','EDA']]


    # Make sure data has a sample rate of 8Hz

    data = interpolateDataTo8Hz(data,sampleRate,startTime)


    # Get the filtered data using a low-pass butterworth filter (cutoff:1hz, fs:8hz, order:6)

    data['filtered_eda'] = butter_lowpass_filter(data['EDA'], 1.0, 8, 6)


    return data

```



```

def interpolateDataTo8Hz(data,sample_rate,startTime):
    if sample_rate<8:
        # Upsample by linear interpolation
        if sample_rate==2:
            data.index = pd.date_range(start=startTime, periods=len(data), freq='500L')
        elif sample_rate==4:
            data.index = pd.date_range(start=startTime, periods=len(data), freq='250L')
            data = data.resample("125L").mean()
    else:
        if sample_rate>8:
            # Downsample
            idx_range = list(range(0,len(data))) # TODO: double check this one
            data = data.iloc[idx_range[0::int(int(sample_rate)/8)]]
            # Set the index to be 8Hz
            data.index = pd.date_range(start=startTime, periods=len(data), freq='125L')

        # Interpolate all empty values
        data = interpolateEmptyValues(data)
    return data

```

```

def interpolateEmptyValues(data):
    cols = data.columns.values
    for c in cols:
        data.loc[:, c] = data[c].interpolate()

    return data

```

```

def butter_lowpass(cutoff, fs, order=5):

```

```
# Filtering Helper functions
```

```
nyq = 0.5 * fs
```

```
normal_cutoff = cutoff / nyq
```

```
b, a = scipy.signal.butter(order, normal_cutoff, btype='low', analog=False)
```

```
return b, a
```

```
def butter_lowpass_filter(data, cutoff, fs, order=5):
```

```
    # Filtering Helper functions
```

```
    b, a = butter_lowpass(cutoff, fs, order=order)
```

```
    y = scipy.signal.lfilter(b, a, data)
```

```
    return y
```